

Report no. : VI-R2020239S-1
Date : 6 November 2020

Sample Number : VI-2020239S-1
Client Name : UNSW Australia
Client Address : UNSW Sydney, NSW 205, Australia
Sample Receiving Date : 27 October 2020
Testing Period : 27 October 2020 – 6 November 2020
Sample Name : FSI-Defend

General Comment:


Sample was received in ambient condition. As agreed by client, 100 mg/L of sample in dilution water was used as testing sample, and the measured testing sample temperature was 21.5°C.

Test Requested

Determination of the acute toxicity to zebrafish eggs (*Danio rerio*) ISO 15088:2007

Test Results : Please refer to next page

Approved signatory:



Dr Xueping Chen

Chief Technology Officer, UK & EU Registered Toxicologist, Chartered Biologist

Hong Kong Accreditation Service (HKAS) has accredited this laboratory (Reg. No. – TEST) under the Hong Kong Laboratory Accreditation Scheme (HOKLAS) for specific laboratory activities as listed in the HOKLAS directory of accredited laboratories.

The copyright of this report is owned by Vitargent (International) Biotechnology Limited. This report shall only be reproduced in full unless prior written approval is obtained from the Head of the Vitargent laboratory.

This Test Report is issued by the Company subject to its General Conditions of Service printed overleaf or available on request. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues defined therein. Unless otherwise stated the results shown in this test report refer only to the samples tested. This test report cannot be reproduced, except in full, without prior written permission of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this report is unlawful and offenders may be prosecuted to the fullest extent of the law.

VITARGENT (INTERNATIONAL) BIOTECHNOLOGY LIMITED

Unit 516, 5/F., Biotech Centre 2 (Bldg 12), No. 11 Science Park West Avenue, Hong Kong Science Park, Shatin, HONG KONG

香港新界沙田香港科技大道西 11 號生物科技中心 2 座 5 樓 516 室

Tel 電話: (852) 2393-2343

Fax 傳真: (852)2393-7922

www.vitargent.com

Report no. : VI-R2020239S-1
Date : 6 November 2020

Test Results

<u>Sample No.</u>	<u>pH</u>	<u>Oxygen concentration</u>	<u>LIDegg</u>
VI-2020239S-1	7.08*	95.9%	1

*Remarks: pH of sample water was adjusted before test according to the method ISO 15088:2007.

LIDegg: Lowest Ineffective Dilution within a test batch in which at least 90% of the fish eggs do not show any effect according to international standard – ISO 15088:2007.

Disclaimer: The test results related only to the submitted test item with the sample number and receiving date specified on the front page.

**** End of Test Report****

Hong Kong Accreditation Service (HKAS) has accredited this laboratory (Reg. No. – TEST) under the Hong Kong Laboratory Accreditation Scheme (HOKLAS) for specific laboratory activities as listed in the HOKLAS directory of accredited laboratories.

The copyright of this report is owned by Vitargent (International) Biotechnology Limited. This report shall only be reproduced in full unless prior written approval is obtained from the Head of the Vitargent laboratory.

This Test Report is issued by the Company subject to its General Conditions of Service printed overleaf or available on request. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues defined therein. Unless otherwise stated the results shown in this test report refer only to the samples tested. This test report cannot be reproduced, except in full, without prior written permission of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this report is unlawful and offenders may be prosecuted to the fullest extent of the law.

VITARGENT (INTERNATIONAL) BIOTECHNOLOGY LIMITED

Unit 516, 5/F., Biotech Centre 2 (Bldg 12), No. 11 Science Park West Avenue, Hong Kong Science Park, Shatin, HONG KONG

香港新界沙田香港科技大道西 11 號生物科技中心 2 座 5 樓 516 室

Tel 電話: (852) 2393-2343

Fax 傳真: (852)2393-7922

www.vitargent.com